

S-Type Optical

The New Platform for Optical Quantum Technologies
and Material Science



S-Type Optical

Compact Sub-Kelvin Cryostat with Free-beam Access

Optical access

Large free-beam access down to the mK stage in a vertical direction.

Proprietary cryogenic interface

Allowing simple integration of experiments prepared outside the cryostat.

Continuous and cryogen free operation

Continuous operation at 300 mK independent of helium-3 supply.

Compact platform

The cryostat and all control equipment can be mounted in a 19" rack.



Continuous Cooling
at **300 mK**



Free-beam Access
Ø > 20 mm





Vibration Level
< 200 nm



Sample Magnet
3 T || optical axis


Cryogen-free continuous operation through cADR

In its standard configuration, the S-Type Optical uses two ADR units to generate continuous magnetic cooling (cADR). The top figure shows the temperature of the ADR units running in continuous mode at 300 mK. While the assisting ADR unit cycles between the 4 K main heat bath (provided by the cryocooler) and a temperature below the target temperature, the second unit controls the sample temperature. As a result, the sample stage can maintain a constant temperature of 300 mK with a typical temperature stability of less than 0.1%, and slightly reduced stability of 4% when activating its heat switch (HS) to initialize the regeneration. The bottom figure shows the temperature control, stepping the sample stage temperature between 200 mK and 1.5 K.



Specs					
System size (w x l x h)	cryostat & electronics	94 x 94 x 208 cm	Temperature stability	typical	< 0.1%
	compressor	54 x 45 x 50 cm		while switching HS	< 4%
Weight	cryostat & electronics	< 600 kg	Cooling power	@300 mK @1 K	15 μ W 250 μ W
Cooldown time	cryostat	32 hrs	Free-beam access	diameter	> 20 mm
Available sample space	diameter	47 mm	Vibration	@1500 Hz	< 200 nm
	height	100 mm		switching HS	< 1.5 μ m

See all specs ►
kiutra.com/s-type-optical



Options

► **3 T sample magnet** smooth bipolar operation to study magnetic properties

► **Wiring**
40 DC and up to 4 RF to the sample stage